

MULTIDEX® GEL/POWDER MALTODEXTRIN

	DO NOT RESTERILIZE
	SINGLE PATIENT USE
STERILE R	STERILIZED WITH IRRADIATION
	NOT MADE WITH NATURAL RUBBER LATEX
RX ONLY	FEDERAL U.S.A. LAW RESTRICTS THIS DEVICE TO SALE OR USE BY OR ON THE ORDER OF A PHYSICIAN OR PROPERLY LICENSED PRACTITIONER.

IMPORTANT INFORMATION

Please read all instructions, warnings, and precautions before use. Correct application is essential for proper functioning of the product. Use only on the person it was provided to by a healthcare professional and only for the use it was intended.

INTENDED USE

Multidex® wound filler is intended to be used on infected and non-infected wounds.

INDICATIONS FOR USE

- Multidex® Gel and Powder wound fillers are indicated for the following wound types:
 - All wounds both infected and non infected
 - Dermal ulcers - partial and full thickness wounds
 - Leg ulcers
 - Pressure ulcers
 - Other exudative lesions
 - Diabetic ulcers
 - Abdominal wounds
 - Superficial wounds
 - Lacerations, cuts, and abrasions
 - Donor sites
 - Second degree burns
- Multidex® Powder wound filler is indicated for moist wounds.
- Multidex® Gel wound filler is indicated for moist to dry wounds.

CONTRAINDICATIONS

- Multidex® Gel and Powder wound fillers are contraindicated for use on 3rd degree burns.

WARNINGS

- Do not use if package is damaged or opened; the seal under the cap is broken or missing; or if the gel is dry
- Do not use if the patient is allergic to this product or any of its ingredients. If the patient has an allergic reaction, cleanse area with sterile saline or warm water, and

discontinue use.

- Do not use if the product is expired. Product may be used for up to 30 days after opening or until gel becomes dry, whichever comes first. Do not contaminate the product applicator with foreign particulates, such as dirt or dust.
- Do not share among multiple users.

CAUTIONS

- Do not use in or near the eye.
- Do not inhale.
- Multidex® wound filler may cause a slight stinging which should stop soon after application. If stinging persists irrigate wound with sterile saline and discontinue use.
- Inappropriate use or use of an inappropriate secondary dressing may lead to adhesion of wound dressing to the tissue. If adhesion occurs, moisten dressing with sterile saline or warm water to dissolve Multidex® wound filler.

DIRECTIONS FOR USE

PREPARATION OF SITE

- Necrotic tissue should be debrided according to clinician or hospital protocol.
- The site should be irrigated liberally with a sterile physiological 0.9% normal saline or a balanced salt solution.

APPLICATION OF MULTIDEX® WOUND FILLER POWDER OR GEL

- After irrigation, apply the Multidex® wound filler.
 - FOR SHALLOW WOUNDS:** layer about 1/4" (6mm) thick over entire wound site.
 - FOR DEEPS WOUNDS:** fill wound site to surface, taking care to fill all undermined areas.
- Cover with a non-adherent dressing, such as MultiPad™, Sofsorb®, Covaderm Plus® or Polyderm Border™. If necessary, tape in place or use Fluftex™ roll gauze or Stretch Net™ to secure dressing.

REMOVAL OF MULTIDEX® WOUND FILLER

- Remove non-adherent dressing with care. If dressing adheres to wound, soak with saline for several minutes before removing so as not to disturb fragile healing tissue.
- Flush site liberally, but gently, with a sterile, physiological irrigating solution to remove debris. This will leave newly formed granulation tissue undisturbed.
- Removal of all Multidex® wound filler is not necessary to continue treatment. Some of the Multidex® wound filler will remain and will blend with more Multidex® wound filler on subsequent applications.
- Follow hospital protocol for disposal of contaminated wound dressings.

Product labels found on boxes may contain a suffix following the product number (Ex. 46-701-**EU**.) This information is for Internal Use Only and does not impact the product. The product itself does not change. Please refer to table below for suffix definitions.

-EU	EUROPEAN SKU
G	SINGLE TUBE
-1	SINGLE BOX

STORAGE AND TRANSPORT CONDITIONS

	KEEP DRY
	KEEP AWAY FROM SUNLIGHT
	DO NOT USE IF PACKAGE IS DAMAGED

WARRANTY

DeRoyal® products are warranted for one hundred twenty (120) days from the date of shipment from DeRoyal as to product quality and workmanship. **DEROYAL'S WRITTEN WARRANTIES ARE GIVEN IN LIEU OF ANY IMPLIED WARRANTIES, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.**